INNOLUX

RED MASK

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asiakaspalvelu@innolux.fi

www.innolux.fi

RED, BLUE, YELLOW AND NEAR-INFRARED LIGHT

Although there is extensive research on the effects of different light wavelengths, this light mask is not a medical device. It is not intended for the diagnosis, treatment, or prevention of diseases, and general research findings cannot be directly applied to this specific device. The product is designed to support general well-being and skincare. In unclear situations, we recommend consulting a healthcare professional.

Wavelength determines how light interacts with the skin and tissues. Red light and near-infrared light penetrate deeper into the skin layers and can reach the mitochondria of cells, which are responsible for cellular respiration. Cellular energy production (ATP) can be activated when cells receive more nutrients and their function becomes more efficient. Yellow and blue light primarily affect the surface layers of the skin.

RED LIGHT (630 nm)

Penetrates the surface layers of the skin and can interact with the mitochondria of cells, which may activate cellular energy production (ATP) and enhance cellular function (photobiomodulation).

- May contribute to skin renewal and improve its appearance
- May support local blood circulation
- May promote recovery of the skin and tissues

NEAR-INFRARED LIGHT (850 nm)

Invisible to the eye. Penetrates deeper into the skin layers than other light spectra and can reach muscle and joint areas. Near-infrared light can interact with the mitochondria of cells and activate cellular energy production (photobiomodulation).

- May promote tissue recovery
- May support local blood circulation
- May contribute to relaxation of the treated area

YELLOW / ORANGE LIGHT (590 nm)

Is mainly absorbed in the surface layers of the skin and does not penetrate tissues as deeply as red or near-infrared light. May affect the function of surface skin cells and microcirculation.

- May have a calming effect on the skin
- May help balance skin redness and tone variations
- May support lymphatic circulation

BLUE LIGHT (460 nm)

Primarily affects the surface layers of the skin. Blue light may have antibacterial properties.

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MAINTENANCE AND WARRANTY

Keep the device clean from dust and dirt. Wipe the mask with a slightly damp cloth or an antibacterial cleaning wipe (note: make sure the mask is not connected to the controller).

Important:

- Keep the device out of reach of children.
- Do not use abrasive or corrosive substances for cleaning.
- Do not use the device in sauna areas or near heat-emitting appliances.
- Do not use the device in damp environments or near water.
- Do not use a damaged device.
- Do not attempt to repair a damaged device yourself.

The LED light sources are not replaceable.

The entire light fixture must be replaced once the light sources reach the end of their lifespan.

Warranty: 2 years

FREQUENTLY ASKED QUESTIONS

HOW SHOULD THE DEVICE BE USED?

We recommend a usage time of 10–20 minutes per session, 2–3 times per week. Use on consecutive days is not recommended. Excessive use is not dangerous, but may reduce results. The best results are achieved by adhering to the recommended usage times and amounts.

CAN THE DEVICE BE USED WITH SKINCARE PRODUCTS?

We recommend using the device on clean skin.

WHEN DO THE EFFECTS OF USING THE DEVICE BECOME VISIBLE?

The first changes may become visible within 2–4 weeks. The first visible change may be brighter and more even-looking skin. More noticeable results (softening of fine lines, improved firmness) may appear after 6–8 weeks. Any potential increase in collagen is gradual, so maximum benefit may develop over time. Results require regular and proper use of the device.

CONTACT INFORMATION

If you have any questions about new Innolux red light equipment, please contact our Customer Service by email or phone. We will be happy to help you! More information is also available on our website.

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